

Ebulliometric Determination of %Alcohol in Wine

Ebulliometry involves the measurement of the decrease in the boiling point (BP) caused by the presence of alcohol in a wine sample. Because the boiling point of a solution also varies with atmospheric pressure, the effects of ethanol are evaluated relative to the boiling point of a control sample of distilled water (determined under identical conditions). The relationship between Δt ($BP_{\text{water}} - BP_{\text{wine}}$) and the %alcohol in the original wine sample is depicted in Figure E-1.

Procedure

Step 1: Determination of Boiling Point of Wine Sample (1:1 dilution)

- a. Add 100ml of distilled water to boiling flask
- b. Add 100ml of wine sample to boiling flask (1:1 dilution)
- c. Insert temperature probe into central port of boiling flask (see Fig E-2).
Tip of temperature probe should be submerged approximately 1 inch.
- d. Attach condenser to boiling flask.
Connect hoses and turn-on cooling water (slow but steady flow rate).
- e. Turn-on heating mantle to highest setting.
- f. Monitor temperature increases.
Begin recording temperature at 15-sec intervals once it exceeds 90°C.
- g. Boiling point is reached when temperature remains unchanged for 3 consecutive samples.
- h. Record boiling point for diluted wine sample (BP_{wine}).

Step 2: Determination of Boiling Point for Distilled Water (Atmospheric Control)

- a. Allow apparatus to cool.
- b. Thoroughly rinse out contents of the boiling flask.
- c. Add 200ml of distilled water to the boiling flask.
- d. Determine boiling point using same procedure described in Step 1 above. Record BP_{water} .
Note: Cooling condenser not needed for Step 2 as there is no risk of alcohol evaporation.

Step 3. Compute %Alcohol Estimate

- a. Determine difference in boiling points (Δt) as follows:
$$\Delta t = BP_{\text{water}} - BP_{\text{wine}}$$
- b. Use Equation E-1 and Δt to determine %Alcohol in the original wine sample.
Note: Figure/Equation E-1 has been adjusted to account for the 1:1 dilution of the wine.

Equation E-1.

$$\% \text{Alcohol} = 0.435 + 1.6687(\Delta t) + 0.1234(\Delta t^2)$$

Note:

Residual sugar levels greater than 4% in the original wine sample may interfere with estimates of %alcohol. The effects of high residual sugar levels can be corrected by applying the following formula:

$$\% \text{alcohol} = \text{est. \%alcohol} - (\% \text{residual sugar} * 0.025)$$




Figure E-1.

Estimated percent alcohol using a 1:1 diluted sample of wine as a function of the difference in the boiling points of the sample versus a distilled water control.

Equation E-1.

$$\% \text{Ethanol} = 0.435 + 1.6687(\Delta t) + 0.1234(\Delta t^2)$$

Ebulliometry Equipment List

100 ml pipette or graduated cylinder for collecting proper volume of wine sample

Electromantle 500-ml heating mantle

500-ml 2-neck boiling flask (24/40)

Leibig condenser (24/40)

IMC 2400 precision digital thermometer

IMC 8011-2009 temperature probe

Drilled #6 rubber stopper (for temperature probe)

Laboratory stand and miscellaneous clamps

Hoses for condenser cooling water




Figure E-2. Ebulliometry Apparatus